

Data sheet Dilution buffer B (10x)

Cat. No.:	CT364
Specificity:	To be used for preparing the working solution of detector antibodies and GABA conjugate in the ELISPOT assay (performed on PVDF membrane-bottomed plates).
Application:	Prior to use, the solution should be diluted 10 times with PBS. Once prepared, the buffer is ready to dilute detector antibodies and GABA conjugate.
Contents:	8 ml. Preservative free.
Sterility:	Membrane filtered (0.2 µm).
Packaging:	10 ml clear glass bottle.
Storage:	Refrigerated. The solution can be safely stored at 4°C for a defined length of time (expiry date indicated on vial). After opening, the product expires after 6 months when kept sterile.
General:	A clear light-yellow solution. Do not use if the solution is or has become turbid.
Safety:	The product is not classified as hazardous according to Regulation (EC) no. 1272/2008 and to Directive 1999/45/EC and its amendments. Consult the material safety data sheet (MSDS) for further information regarding hazards and safe handling practices. The MSDS can be found on our website www.ucytech.com .

