

## Data sheet Cytokine Stabilization buffer (CSB)

<b>Cat. No.:</b>	<b>CT373</b>
<b>Specificity:</b>	Inhibits cytokine degradation in plasma and serum samples of human and animal origin during incubation at 37°C.
<b>Application:</b>	Add to non-diluted serum or plasma samples at 1/20 volume just prior to cytokine analysis, using U-CyTech ELISA systems.
<b>Contents:</b>	5 ml. Supplied in PBS supplemented with 1% bovine serum albumin (BSA). Preservative free.
<b>Sterility:</b>	Membrane filtered (0.2 µm).
<b>Packaging:</b>	10 ml clear glass bottle.
<b>Storage:</b>	Refrigerated. The solution can be safely stored at 4°C for a defined length of time (expiry date indicated on vial). After opening, the product expires after 6 months when kept sterile.
<b>General:</b>	A clear light-yellow solution. Do not use if the solution has become turbid.
<b>Safety:</b>	The product is not classified as hazardous according to Regulation (EC) no. 1272/2008 and to Directive 1999/45/EC and its amendments. Consult the material safety data sheet (MSDS) for further information regarding hazards and safe handling practices. The MSDS can be found on our website <a href="http://www.ucytech.com">www.ucytech.com</a> .

