

**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Monkey IL-17A ELISPOT kit

Catalogue no.: CT401-PR2 (2-plate format)  
CT401-PR5 (5-plate format)

### Validated for:

Rhesus macaque and Baboon



Stimulus: none ( $2 \times 10^5$  cells/well)



Stimulus: Concanavalin A ( $5 \times 10^4$  cells/well)

Example of IL-17A specific spots produced by rhesus macaque PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY   
ANIARA DIAGNOSTICA LLC  
+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com