

U-CyTech BV


Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-17F ELISPOT kit


Catalogue no.: CT418-T2 (2-plate format)
CT418-T5 (5-plate format)

Validated for:

Human



Stimulus: none (2×10^5 cells/well)



Stimulus: PMA/ionomycin (1×10^5 cells/well)

Example of IL-17F specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 16h.

A preincubation step (indirect procedure) is recommended for spot development.



PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com