


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Mouse IL-6 ELISPOT kit


Catalogue no.: CT436-T2 (2-plate format)
CT436-T5 (5-plate format)

Validated for:

Mouse



Stimulus: none (1×10^5 cells/well)



Stimulus: Concanavalin A (1×10^5 cells/well)

Example of IL-6 specific spots produced by mouse (CD-1) splenocytes.
Preincubation: 24h, incubation ELISPOT: 20h.



PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com