


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-17A ELISA kit

Catalogue no.: CT516A (5-plate format)

Validated for:

Human



Example of a representative Human IL-17A standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.
Sensitivity: 2 pg/ml.
Specificity: natural and recombinant human IL-17A.
Calibration: recombinant human IL-17A.
Type of sample: serum, plasma and culture supernatant.

References list:

Beigier-Bompadre, M., Moos, V., Belogolova, E., Allers, K., Schneider, T., Churin, Y., Ignatius, R., Meyer, T.F. and Aebischer, T. 2011. J. Infect. Dis. 204: 1339-1348

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 140820