

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Human IL-23 ELISA kit

Catalogue no.: CT517A (5-plate format)

Validated for:

Human



Example of a representative Human IL-23 standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.

Sensitivity: 2 pg/ml.

Specificity: natural and recombinant human IL-23.

Calibration: recombinant human IL-23.

Type of sample: culture supernatant.

No cross-reactivity or interference was observed with IL-12p40 or IL-12p70.

References list:

Banik, S., Renner Viveros, P., Seeber, F., Klotz, C., Ignatius, R. and Aebischer, T. 2013. Infect. Immun. 81: 2309-2317

Fehlings, M., Drobbe, L., Moos, V., Renner Viveros, P., Hagen, J., Beigier-Bompadre, M., Pang, E., Belogolova, E., Churin, Y., Schneider, T., Meyer, T.F., Aebischer, T. and Ignatius, R. 2012. Infect. Immun. 80: 2724-2734

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 140820