


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-17F ELISA kit

Catalogue no.: CT518A (5-plate format)

Validated for:

Human



Example of a representative Human IL-17F standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.
Sensitivity: 4 pg/ml.
Specificity: natural and recombinant human IL-17F.
Calibration: recombinant human IL-17F.
Type of sample: serum, plasma and culture supernatant.

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 140820