

U-CyTech BV

Yalelaan 48
3584 CM Utrecht
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

Typical data Human IL-31 ELISA kit

Catalogue no.: CT520A (5-plate format)

Validated for:

Human



Example of a representative Human IL-31 standard curve using a 4-parameter logistic fit.

Recommended assay range: 32-2048 pg/ml.

Sensitivity: 32 pg/ml.

Specificity: natural and recombinant human IL-31.

Calibration: recombinant human IL-31.

Type of sample: culture supernatant.

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797

7768 Service Center Drive
West Chester, OH 45069

info@aniara.com

Version 160708