


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-27 ELISA kit

Catalogue no.: CT524A (5-plate format)

Validated for:

Human



Example of a representative Human IL-27 standard curve using a 4-parameter logistic fit.

Recommended assay range: 5-320 pg/ml.

Sensitivity: 5 pg/ml.

Specificity: natural and recombinant human IL-27.

Calibration: recombinant human IL-27.