

U-CyTech BV

Yalelaan 48  
3584 CM Utrecht  
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

## Typical data Human IL-1 $\beta$ ELISA kit

Catalogue no.: CT526A (5-plate format)

### Validated for:

Human



Example of a representative Human IL-1 $\beta$  standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.

Sensitivity: 2 pg/ml.

Specificity: natural and recombinant human IL-1 $\beta$ .

Calibration: recombinant human IL-1 $\beta$ .

Type of sample: serum, plasma and culture supernatant.

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com

Version 160216