

Data sheet Human IL-12p70 ELISPOT antibody pair; 20-plate format

Cat. No.: CT650-20

Coating antibodies (1 vial)

Product: Monoclonal antibody to human interleukin 12p70 (IL-12p70)
(does not bind the p40 or p35 subunits)

Isotype: Mouse IgG₁

Production: *In vitro* using serum free medium

Purification: DEAE ion exchange chromatography

Contents: Each vial contains sufficient material for coating of twenty 96-well ELISPOT plates

Buffer: Prior to lyophilization: 1.0 ml PBS + 125 mM trehalose

Application: Coating antibody in an ELISPOT system

Reconstitution: Dissolve the contents of the vial by injection of 1.0 ml distilled water into the vial and dilute 100 times in PBS. The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 50 µl/well).

Detection antibodies (1 vial)

Product: Biotinylated polyclonal antibody to human interleukin 12p70 (IL-12p70)

Isotype: Rabbit IgG

Purification: Ammonium sulphate precipitation and protein A-affinity chromatography

Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)

Contents: Each vial contains sufficient material for twenty 96-well ELISPOT plates

Buffer: Prior to lyophilization: 2.0 ml PBS + 1% BSA + 125 mM trehalose

Application: Detection antibody in an ELISPOT system

Reconstitution: Dissolve the contents of the vial by injection of 2.0 ml distilled water into the vial and dilute 100 times in Dilution buffer (see Technical Data Sheet). The total amount of one vial is sufficient for twenty 96-well ELISPOT plates (1920 determinations; 100 µl/well).

General

Specificity: Validated for detecting human IL-12

Sterility: Membrane filtered (0.2 µm)

Stability: The lyophilized products are stable for at least one year at 4°C (expiry date is indicated on the vials).
After reconstitution, the antibodies are stable for several months at 4°C (if kept sterile) or for minimal one year at -20°C.

For research use only

PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 111025