

U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Data sheet

Monoclonal antibody to human/mouse interleukin 5 (IL-5)

Cat. No.:	CT694
Clone:	TRFK5
Isotype:	Rat IgG ₁
Production:	<i>In vitro</i> using serum free medium.
Purification:	Ion exchange chromatography.
Packaging:	Lyophilized and vacuum-packed.
Contents:	0.5 mg/vial
Buffer:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Specificity:	Binds with high efficiency to natural and recombinant human and mouse IL-5.
Sterility:	Membrane filtered (0.2 µm).
Reconstitution:	Dissolve the contents of the vial by injection of 0.5 ml sterile distilled water.
Stability:	Lyophilized product is stable for at least one year at 4°C (expiry date is indicated on the vial). After reconstitution, the contents can be safely stored at 4°C for one month or for one year at -20°C. Add 0.02% sodiumazide to prevent bacterial growth.
Quantitation:	Antibody concentration was determined by absorbance, taking A ₂₈₀ =1.4 for a 1 mg/ml solution.
Application:	Coating antibody in ELISA Intracellular staining <i>In vitro</i> neutralization Western blot analysis



For research use only

PROUDLY DISTRIBUTED BY 
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com

Version 150706