

## Data sheet **Mouse IL-2 ELISA antibody pair; 10-plate format**

Cat. No.: CT762-10

### Coating antibodies (2 vials)

**Product:** Monoclonal antibody to mouse interleukin 2 (IL-2)  
**Isotype:** Rat IgG<sub>2a</sub>  
**Production:** *In vitro* using Hybridoma-SFM medium  
**Purification:** Protein G chromatography  
**Buffer:** Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose  
**Application:** Inject 0.25 ml distilled water into the vial and dilute 100 times in PBS. The content of one vial is sufficient for five 96-well ELISA plates (480 determinations; 80 µl/well).

### Detection antibodies (2 vials)

**Product:** Biotinylated monoclonal antibody to mouse interleukin 2 (IL-2)  
**Isotype:** Rat IgG<sub>2b</sub>  
**Production:** *In vitro* using serum free medium  
**Purification:** Protein G chromatography  
**Labeling:** With Biotin-7-NHS (N-hydroxysuccinimide)  
**Buffer:** Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose  
**Application:** Inject 0.5 ml distilled water into the vial and dilute 100 times in PBS + 0.5% BSA + 0.05% Tween-20. The content of one vial is sufficient for five 96-well ELISA plates (480 determinations; 100 µl/well).

### Standards (5 vials)

**Product:** Recombinant mouse interleukin 2 (IL-2)  
**Application:** Cytokine standard in an ELISA system  
**Reconstitution:** Dissolve the contents of one vial by injection of 0.5 ml distilled water into the vial. Use immediately.

### Conjugate (2 vials)

**Product:** SPP conjugate (Streptavidin-HRP)  
**Application:** Inject 0.5 ml distilled water into the vial and dilute 100 times in PBS + 0.5% BSA + 0.05% Tween-20. The content of one vial is sufficient for five 96-well ELISA plates (480 determinations; 100 µl/well).  
The product should be used in combination with TMB substrate.

### General

**Sensitivity:** 4 pg/ml  
**Specificity:** Validated for detecting natural and recombinant mouse IL-2  
**Sterility:** Membrane filtered (0.2 µm)  
**Stability:** Lyophilized SPP conjugate is stable for at least one year at -20°C in the dark, the other lyophilized products are stable for at least one year at 4°C. After reconstitution, the antibodies are stable for one year at 4°C (if kept sterile) and SPP for minimal one year at -20°C in the dark. The reconstituted standard preparation should be used immediately.

For research use only