

U-CyTech BV

Yalelaan 48  
3584 CM Utrecht  
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

## Typical data Mouse IgG<sub>2a</sub> B cell ELISPOT kit


Catalogue no.: CT792-T5 (5-plate format)

### Validated for:

Mouse



Total number of IgG<sub>2a</sub> secreting cells  
1x10<sup>5</sup> cells/well (Assay I)  
Stimulus preincubation: none



Total number of IgG<sub>2a</sub> secreting cells  
2x10<sup>4</sup> cells/well (Assay I)  
Stimuli preincubation: IL-2 + R848

Example of IgG<sub>2a</sub> specific spots produced by mouse (Balb/c) splenocytes.

B cell activation: 5x10<sup>6</sup> cells/ml, 2 days preincubation without stimulus (left) or in the presence of recombinant IL-2 and R848 (right).

Incubation ELISPOT: 6 hours.

PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com

Version 140530