

**U-CyTech BV**  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com


## Typical data Human IFN- $\gamma$ /IL-17A FluoroSpot kit

Catalogue no.: CT855-F2 (2-plate format)

### Validated for:

Human

Example of IFN- $\gamma$ /IL-17A specific spots produced by  $1 \times 10^4$  human PBMC:



In green: IFN- $\gamma$ , in red: IL-17A, in yellow: IFN- $\gamma$ /IL-17A.

Preincubation: 42h, incubation FluoroSpot: 20h.  
Stimulation: PMA/ionomycin.



PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com