

U-CyTech BV  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com


## Typical data Monkey IFN- $\gamma$ /TNF- $\alpha$ FluoroSpot kit

Catalogue no.: CT901-F2 (2-plate format)

### Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey

Example of IFN- $\gamma$ /TNF- $\alpha$  specific spots produced by  $2.5 \times 10^4$  rhesus macaque PBMC:



In green: IFN- $\gamma$ , in red: TNF- $\alpha$ , in yellow: IFN- $\gamma$ /TNF- $\alpha$ .

Preincubation: 42h, incubation FluoroSpot: 20h.  
Stimulation: PMA/ionomycin.



PROUDLY DISTRIBUTED BY

**ANIARA**  
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797  
7768 Service Center Drive  
West Chester, OH 45069  
info@aniara.com