

U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com


Typical data Monkey IFN- γ /IL-2 FluoroSpot kit

Catalogue no.: CT902-F2 (2-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey

Example of IFN- γ /IL-2 specific spots produced by 2.5×10^4 rhesus macaque PBMC:



In green: IFN- γ , in red: IL-2, in yellow: IFN- γ /IL-2.

Preincubation: 42h, incubation FluoroSpot: 20h.
Stimulation: PMA/ionomycin.



PROUDLY DISTRIBUTED BY

ANIARA
ANIARA DIAGNOSTICA LLC

+1 (513) 770-1991 | +1 (866) 783-3797
7768 Service Center Drive
West Chester, OH 45069
info@aniara.com