


Typical data Monkey IFN- γ ELISPOT kit


Catalogue no.: CT121-PB2 (2-plate format)
CT121-PB5 (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey



Stimulus: none (4×10^5 cells/well)



Stimulus: Concanavalin A (1×10^4 cells/well)

Example of IFN- γ specific spots produced by rhesus macaque PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

