


## Typical data Monkey IFN- $\gamma$ ELISPOT kit

Catalogue no.: CT121-PR2 (2-plate format)  
CT121-PR5 (5-plate format)

### Validated for:

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey



Stimulus: none ( $4 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $1 \times 10^4$  cells/well)

Example of IFN- $\gamma$  specific spots produced by rhesus macaque PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.

### References list:

Zhao, L., Liu, B., Ren, J., Feng, J., Pang, Z., Gao, J., Zhang, H., Tan, W., Tian, H. and Ruan L. 2011. *Virology* 15: 302