


## Typical data **Monkey IFN- $\gamma$ ELISPOT kit**


**Catalogue no.:** CT121-T2 (2-plate format)  
CT121-T5 (5-plate format)

### **Validated for:**

Rhesus macaque, Cynomolgus monkey, Pig-tailed macaque, Baboon, African green monkey



Stimulus: none ( $5 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $2.5 \times 10^4$  cells/well)

Example of IFN- $\gamma$  specific spots produced by rhesus macaque PBMC.  
Preincubation: 42h, incubation ELISPOT: 5h.

