


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com


Typical data **Monkey Perforin ELISPOT kit**

Catalogue no.: CT136-PR2 (2-plate format)
CT136-PR5 (5-plate format)

Validated for:
Rhesus macaque, Barbary macaque, Baboon



Stimulus: none (2×10^5 cells/well)



Stimulus: PMA/ionomycin (1×10^5 cells/well)

Example of perforin specific spots produced by baboon PBMCs.
Preincubation: 48h, incubation ELISPOT: 20h.

