


Typical data Monkey IL-8 ELISA kit

Catalogue no.: CT151A (5-plate format)

Validated for:

Rhesus macaque, Barbary macaque, Baboon



Example of a representative Monkey IL-8 standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.

Sensitivity: 8 pg/ml.

Calibration: natural rhesus macaque IL-8.

Type of sample: plasma, serum and culture supernatant.

References list:

Jacquelin, B., Petitjean, G., Kunkel, D., Liovat, A.-S., Jochems, S.P., Rogers, K.A., Ploquin, M.J., Madec, Y., Barré-Sinoussi, F., Dereuddre-Bosquet, N., Lebon, P., Le Grand, R., Villinger, F. and Müller-Trutwin, M. 2014. PLoS Pathog. 10: e1004241