

**U-CyTech BV**

Yalelaan 48  
3584 CM Utrecht  
The Netherlands

P +31.30.253 5960

F +31.30.253 9344

INFO@ucytech.com


www.ucytech.com

## Typical data Monkey IP-10 ELISA kit

Catalogue no.: CT157A (5-plate format)

### Validated for:

Rhesus macaque



Example of a representative Monkey IP-10 standard curve using a 4-parameter logistic fit.

Recommended assay range: 5-320 pg/ml.

Sensitivity: 5 pg/ml.

Calibration: natural rhesus macaque IP-10.

Type of sample: serum, plasma and culture supernatant.