


Typical data Human GM-CSF ELISA kit

Catalogue no.: CT200A (5-plate format)

Validated for:

Human



Example of a representative Human GM-CSF standard curve using a 4-parameter logistic fit.

Recommended assay range: 2-128 pg/ml.
Sensitivity: 2 pg/ml.
Specificity: natural and recombinant human GM-CSF.
Calibration: recombinant human GM-CSF.
Type of sample: serum, plasma and culture supernatant.

References list:

Duan, R., Remeijer, L., van Dun, J.M., Osterhaus, A.D.M.E. and Verjans, G.M.G.M. 2007. Invest. Ophthalmol. Vis. Sci. 48: 277-284