


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human Granzyme B ELISA kit

Catalogue no.: CT211A (5-plate format)

Validated for:

Human



Example of a representative Human Granzyme B standard curve using a 4-parameter logistic fit.

Recommended assay range: 25-1600 pg/ml.

Sensitivity: 25 pg/ml.

Specificity: natural and recombinant human Granzyme B.

Calibration: natural human Granzyme B.

Type of sample: serum, plasma and culture supernatant.