


Typical data Human IFN- γ ELISPOT kit


Catalogue no.: CT230-PB2 (2-plate format)
CT230-PB5 (5-plate format)

Validated for:

Human



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (1×10^4 cells/well)

Example of IFN- γ specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

Chiriva-Interati, M., Mirandola, L., Yu, Y., Jenkins, M.R., Gornati, R., Bernardini, G., Gioia, M, Chiamonte, R., Cannon, M.J., Kast, W.M. and Cobos, E. 2011. *J. Immunother.* 34: 490-499

Correale, P., Botta, C., Cusi, M., Del Vecchio, M., De Santi, M., Gori Savellini, G., Bestoso, E., Apollinari, S., Mannucci, S., Marra, M., Abruzzese, A., Aquino, A., Turriziani, M., Bonmassar, L., Caraglia, M. and Tagliaferri, P. 2012. *Int. J. Cancer* 130: 1577-1589

Ostrov, D.A., Alkanani, A., McDaniel, K.A., Case, S, Baschal, E.E., Pyle, L., Ellis, S., Pöllinger, B., Seidl, K.J., Shah, V.N., Garg, S.K., Atkinson, M.A. Gottlieb, P.A. and Michels, A.W. 2018. *J. Clin. Invest.* 128: 1888-1902

More references can be found in our References Database on our website: www.ucytech.com/references.