


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.30.253 5960
F +31.30.253 9344
INFO@ucytech.com
www.ucytech.com

Typical data Human IL-2 ELISPOT kit


Catalogue no.: CT231-T2 (2-plate format)
CT231-T5 (5-plate format)

Validated for:

Human



Stimulus: none (2×10^5 cells/well)



Stimulus: PMA/ionomycin (2.5×10^4 cells/well)

Example of IL-2 specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 5h.

References list:

- Van Halteren, A.G.S., van Etten, E., de Jong, E.C., Bouillon, R., Roep, B.O. and Mathieu, C. 2002. *Diabetes* 51: 2119-2125
- Kloosterboer, F.M., van Luxemburg-Heijs, S.A.P., Willemze, R. and Falkenburg, J.H.F. 2004. *Hum. Immunol.* 65: 328-339
- Lv, M., Li, Y., Yu, M., Sun, Y., Lin, Z., Qiao, C., Luo, Q., Gu, X., Huang, Y., Feng, J. and Shen, B. 2007. *Int. J. Biochem. Cell Biol.* 39: 1142-1155
- Skowera, A., de Jong, E.C., Schuitemaker, J.H.N., Allen, J.S., Wessely, S.C., Griffiths, G., Kaspenberg, M. and Peakman, M. 2005. *J. Immunol.* 175: 7235-7243