


U-CyTech BV  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Human IL-4 ELISPOT kit


Catalogue no.: CT232-PR2 (2-plate format)  
CT232-PR5 (5-plate format)

### Validated for:

Human



Stimulus: none ( $4 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $2.5 \times 10^4$  cells/well)

Example of IL-4 specific spots produced by human PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.

### References list:

Xue, X-K., Sun, W-J., Hu, J., Yao, H-P. and Cai, Z. 2006. Chin. Med. J. 119: 1038-1041