

U-CyTech BV  
Yalelaan 48  
3584 CM Utrecht  
The Netherlands  
P +31.30.253 5960  
F +31.30.253 9344  
INFO@ucytech.com  
www.ucytech.com

## Typical data Human GM-CSF ELISPOT kit


Catalogue no.: CT241-PR2 (2-plate format)  
CT241-PR5 (5-plate format)

### Validated for:

Human



Stimulus: none ( $4 \times 10^5$  cells/well)



Stimulus: PMA/ionomycin ( $1 \times 10^4$  cells/well)

Example of GM-CSF specific spots produced by human PBMC.  
Preincubation: 42h, incubation ELISPOT: 20h.

