


## Typical data Mouse IL-10 ELISA kit

Catalogue no.: CT307A (5-plate format)

### Validated for:

Mouse



Example of a representative Mouse IL-10 standard curve using a 4-parameter logistic fit.

Recommended assay range: 32-2048 pg/ml.  
 Sensitivity: 32 pg/ml.  
 Specificity: natural and recombinant mouse IL-10.  
 Calibration: recombinant mouse IL-10.  
 Type of sample: serum, plasma and culture supernatant.

### References list:

- Ghotloo, S., Haji Molla Hoseini, M., Alimohammadian, M.H., Khaze, V., Memarnejadian, A. and Rostami, A. 2015. Parasitol. Int. 64: 219-221
- Hong, C., Zhang, T. and Gao, X-M. 2013. Mol. Immunol. 55: 237-246
- Pouw, N., Treffers-Westerlaken E., Kraan J., Wittink F., Ten Hagen T., Verweij J. and Debets R. 2010. Cancer Immunol. Immunother. 59: 921-931
- Willemen, H.L.D.M., Eijkelkamp, N., Garza Carbajal, A., Wang, H., Mack, M., Zijlstra, J., Heijnen, C.J. and Kavelaars, A. 2014. J. Pain 15: 496-506
- Zhao, Y., Liu, M.Y., Wang, X.L., Liu, X.L., Yang, Y., Zou, H.B., Sun, S.M., Yu, L., Rosenthal, B., Shi, H.N., Boireau, P. and Wu, X.P. 2013. Vet. Parasitol. 194: 211-216

More references can be found in our References Database on our website: [www.ucytech.com/references](http://www.ucytech.com/references).