


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Mouse IL-2 ELISA kit

Catalogue no.: CT309A (5-plate format)

Validated for:

Mouse



Example of a representative Mouse IL-2 standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.
Sensitivity: 4 pg/ml.
Specificity: natural and recombinant mouse IL-2.
Calibration: recombinant mouse IL-2.
Type of sample: serum, plasma and culture supernatant.