


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com


Typical data Mouse IL-5 ELISPOT kit

Catalogue no.: CT321-PR2 (2-plate format)
CT321-PR5 (5-plate format)

Validated for:
Mouse



Stimulus: none (2×10^5 cells/well)



Stimulus: PMA/ionomycin (2.5×10^4 cells/well)

Example of IL-5 specific spots produced by mouse (Balb/c) splenocytes.
Preincubation: 42h, incubation ELISPOT: 20h.

