


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human Perforin ELISA kit

Catalogue no.: CT391A (5-plate format)

Validated for:

Human



Example of a representative Human Perforin standard curve using a 4-parameter logistic fit.

Recommended assay range: 64 - 4096 pg/ml.

Sensitivity: 64 pg/ml.

Specificity: natural human Perforin.

Calibration: natural human Perforin.

Type of sample: serum, plasma and culture supernatant.