



Typical data Human IL-17A ELISPOT kit

Catalogue no.: CT416-T2 (2-plate format)
CT416-T5 (5-plate format)

Validated for:
Human



Stimulus: none (4×10^5 cells/well)



Stimulus: PMA/ionomycin (5×10^4 cells/well)

Example of IL-17A specific spots produced by human PBMC.
Preincubation: 42h, incubation ELISPOT: 20h.

References list:

McLaughlin, K.A., Gulati, K., Richardson, C.C., Morgan, D., Bodansky., H.J., Feltbower, R.G. and Christie, M.R. 2014. *J. Immunol.* 193: 4448-4456

van Besouw, N.M., Caliskan, K., Peeters, A.M., Klepper, M., Dieterich, M., Maat, L.P., Weimar, W. Manintveld, O.C. and Baan, C.C. 2015. *J. Heart Lung Transplant.* 34: 933-940