


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Monkey IL-23 ELISA kit

Catalogue no.: CT502A (5-plate format)

Validated for:

Rhesus macaque, Cynomolgus monkey



Example of a representative Monkey IL-23 standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.
Sensitivity: 8 pg/ml.
Calibration: natural rhesus macaque IL-23.
Type of sample: culture supernatant.