


Typical data Human IL-17A ELISA kit

Catalogue no.: CT516A (5-plate format)

Validated for:

Human



Example of a representative Human IL-17A standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.
Sensitivity: 8 pg/ml.
Specificity: natural and recombinant human IL-17A.
Calibration: recombinant human IL-17A.
Type of sample: serum, plasma and culture supernatant.

References list:

Beigier-Bompadre, M., Moos, V., Belogolova, E., Allers, K., Schneider, T., Churin, Y., Ignatius, R., Meyer, T.F. and Aebischer, T. 2011. J. Infect. Dis. 204: 1339-1348