


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human IL-17F ELISA kit

Catalogue no.: CT518A (5-plate format)

Validated for:

Human



Example of a representative Human IL-17F standard curve using a 4-parameter logistic fit.

Recommended assay range: 4-256 pg/ml.

Sensitivity: 4 pg/ml.

Specificity: natural and recombinant human IL-17F.

Calibration: recombinant human IL-17F.

Type of sample: serum, plasma and culture supernatant.