


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human MCP-1 ELISA kit

Catalogue no.: CT521A (5-plate format)

Validated for:

Human



Example of a representative Human MCP-1 standard curve using a 4-parameter logistic fit.

Recommended assay range: 8-512 pg/ml.

Sensitivity: 8 pg/ml.

Specificity: natural and recombinant human MCP-1.

Calibration: recombinant human MCP-1.

Type of sample: culture supernatant.