


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human Angiopoietin-2 ELISA kit

Catalogue no.: CT527A (5-plate format)

Validated for:

Human



Example of a representative Human Angiopoietin-2 standard curve using a 4-parameter logistic fit.

Recommended assay range: 100-6400 pg/ml.

Sensitivity: 100 pg/ml.

Specificity: natural and recombinant human angiopoietin-2.

Calibration: recombinant human angiopoietin-2.

Type of sample: serum, plasma and culture supernatant.