


U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com

Typical data Human IL-21 ELISA kit

Catalogue no.: CT530A (5-plate format)

Validated for:

Human



Example of a representative Human IL-21 standard curve using a 4-parameter logistic fit.

Recommended assay range: 5-320 pg/ml.

Sensitivity: 5 pg/ml.

Specificity: natural and recombinant human IL-21.

Calibration: recombinant human IL-21.

Type of sample: culture supernatant.