

Data sheet **Monkey IL-6 ELISA antibody pair** 10-plate and 20-plate format

Contents

Items	Cat. no. CT714-10 (10-plate format)	Cat. no. CT714-20 (20-plate format)	Storage conditions
Coating antibody	2 vials	4 vials	4 °C
Biotinylated detection antibody	2 vials	4 vials	4 °C
Standards	5 vials	10 vials	4 °C
SPP conjugate	2 vials	4 vials	≤-20 °C in the dark

Please read our ELISA antibody pair Technical Data Sheet for further guidelines and procedures.

Coating antibody

Product: Monoclonal antibody to monkey interleukin 6 (IL-6)
 Isotype: Mouse IgG₁
 Production: *In vitro* using serum free medium
 Purification: Ion exchange chromatography
 Buffer: Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose

Biotinylated detection antibody

Product: Biotinylated polyclonal antibody to monkey interleukin 6 (IL-6)
 Isotype: Goat IgG
 Purification: Protein G-affinity chromatography
 Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)
 Buffer: Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose

Standards

Product: Recombinant interleukin 6 (IL-6)

SPP conjugate

Product: SPP conjugate (Streptavidin-HRP)
 Application: This product should be used in combination with TMB substrate.

General

Sensitivity: 1 pg/ml
 Specificity: Validated for detecting natural rhesus macaque, cynomolgus monkey, pig-tailed macaque, barbary macaque, baboon and African green monkey IL- 6.
 Sterility: Membrane filtered (0.2 µm)
 Content: Lyophilized
 References: Lainée P *et al.* (2005). Crit Care Med 33: 797-805.
 Mehlhop E *et al.* (2002). J Immunol Meth 260: 219-34.
 Wang CY *et al.* (2007). Vaccine 25: 3041-52.

More references can be found at www.ucytech.com/references.