

## Data sheet **Monkey IL-10 ELISA antibody pair** 10-plate and 20-plate format

### Contents

Items	Cat. no. CT715-10 (10-plate format)	Cat. no. CT715-20 (20-plate format)	Storage conditions
Coating antibody	2 vials	4 vials	4 °C
Biotinylated detection antibody	2 vials	4 vials	4 °C
Standards	5 vials	10 vials	4 °C
SPP conjugate	2 vials	4 vials	≤-20 °C in the dark

Please read our ELISA antibody pair Technical Data Sheet for further guidelines and procedures.

### Coating antibody

Product: Monoclonal antibody to monkey interleukin 10 (IL-10)  
 Isotype: Mouse IgG<sub>1</sub>  
 Production: *In vitro* using serum free medium  
 Purification: Ion exchange chromatography  
 Buffer: Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose

### Biotinylated detection antibody

Product: Biotinylated monoclonal antibody to monkey interleukin 10 (IL-10)  
 Isotype: Mouse IgG<sub>1</sub>  
 Production: *In vitro* using serum free medium  
 Purification: Protein A-affinity chromatography  
 Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)  
 Buffer: Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose

### Standards

Product: Recombinant interleukin 10 (IL-10)

### SPP conjugate

Product: SPP conjugate (Streptavidin-HRP)  
 Application: This product should be used in combination with TMB substrate.

### General

Sensitivity: 2 pg/ml  
 Specificity: Validated for detecting natural rhesus macaque, cynomolgus monkey, pig-tailed macaque, baboon and African green monkey IL-10.  
 Sterility: Membrane filtered (0.2 µm)  
 Content: Lyophilized  
 References: Asiedu C *et al.* (2007). Cytokine 40: 183-92.  
 Coe CL *et al.* (2012). Age(Dordr) 34(5): 1169-77.  
 Li M *et al.* (2018). Xenotranspl 25(1).

More references can be found at [www.ucytech.com/references](http://www.ucytech.com/references).