

## Data sheet Human IFN- $\gamma$ ELISA antibody pair 10-plate and 20-plate format

### Contents

Items	Cat. no. CT740-10 (10-plate format)	Cat. no. CT740-20 (20-plate format)	Storage conditions
Coating antibody	2 vials	4 vials	4 °C
Biotinylated detection antibody	2 vials	4 vials	4 °C
Standards	5 vials	10 vials	4 °C
SPP conjugate	2 vials	4 vials	$\leq -20$ °C in the dark

Please read our ELISA antibody pair Technical Data Sheet for further guidelines and procedures.

### Coating antibody

Product: Monoclonal antibody to human interferon gamma (IFN- $\gamma$ )  
 Isotype: Mouse IgG<sub>1</sub>  
 Production: *In vitro* using serum free medium  
 Purification: DEAE ion exchange chromatography  
 Buffer: Prior to lyophilization: 0.25 ml PBS + 125 mM trehalose

### Biotinylated detection antibody

Product: Biotinylated polyclonal antibody to human interferon gamma (IFN- $\gamma$ )  
 Isotype: Rabbit IgG  
 Purification: Ammonium sulphate precipitation, protein A- and ligand-affinity chromatography  
 Labeling: With Biotin-7-NHS (N-hydroxysuccinimide)  
 Buffer: Prior to lyophilization: 0.5 ml PBS + 1% BSA + 125 mM trehalose

### Standards

Product: Recombinant human interferon gamma (IFN- $\gamma$ )

### SPP conjugate

Product: SPP conjugate (Streptavidin-HRP)  
 Application: This product should be used in combination with TMB substrate.

### General

Sensitivity: 2 pg/ml  
 Specificity: Validated for detecting natural and recombinant human IFN- $\gamma$   
 Sterility: Membrane filtered (0.2  $\mu$ m)  
 Content: Lyophilized  
 References: Bobosha K *et al.* (2014). PLoS Negl Trop Dis 8: e2845.  
 Geluk A *et al.* (2012). J Immunol 188: 4782-4791.  
 Iyer AM *et al.* (2007). Hum Pathol 38: 1065-1073.

More references can be found at [www.ucytech.com/references](http://www.ucytech.com/references).