

U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com


Typical data Human IFN- γ /Granzyme B FluoroSpot kit

Catalogue no.: CT849-F2 (2-plate format)

Validated for:

Human

Example of IFN- γ /Granzyme B specific spots produced by 2.5×10^5 human PBMC:



In green: IFN- γ , in red: Granzyme B, in yellow: IFN- γ /Granzyme B.

Preincubation: 42h, incubation FluoroSpot: 20h.
Stimulation: PMA/ionomycin.