

U-CyTech BV
Yalelaan 48
3584 CM Utrecht
The Netherlands
P +31.85.073 1460
info@ucytech.com
www.ucytech.com


Typical data Human IFN- γ /IL-4 FluoroSpot kit

Catalogue no.: CT857-F2 (2-plate format)

Validated for:

Human

Example of IFN- γ /IL-4 specific spots produced by 5×10^4 human PBMC:



In green: IFN- γ , in red: IL-4, in yellow: IFN- γ /IL-4.

Preincubation: 42h, incubation FluoroSpot: 20h.
Stimulation: PMA/ionomycin.

