

Aniara Diagnostica Launches MORE Microsites: 5d-stabilize.com and aprotininproducts.com

Aniara Diagnostica is pleased to announce the addition of more microsites. Easy navigation to your specific product needs! This includes: 5D Protein Stabilizer, 10% Hydrolyzed Porcine Collagen Peptides, and Aprotinin, a multiple inhibitor of many proteases. For Research Use Only. Not for use in diagnostic procedures.

WEST CHESTER, Ohio ([PRWEB](#)) October 27, 2020 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality by launching <https://5d-stabilize.com> and <https://aprotininproducts.com>

This product is for Research Use Only and not for use in Diagnostic procedures.

5D's Protein Stabilizer consists of hydrolyzed pig collagen peptide fraction, from dermal origin, which can be used as a protein stabilizer for lyophilized, liquid formulations or as a blocking agent for coated surfaces. 5D's Protein Stabilizer is free from cartilage, bone and plasma components. The product is sterilized for a long-term storage.

A5D-82411S (50 mL) 5D Protein Stabilizer [click here for more information](#)

A5D-82411A (100 mL) 5D Protein Stabilizer [click here for more information](#)

A5D-82411B (1000mL) 5D Protein Stabilizer [click here for more information](#)

Aprotinin can be used in any application where inhibition of its targeted proteases is needed. Current uses concern protein isolation, removal of reactive proteases through solid phase, inhibition of undesired proteolytic activity in laboratory assays, or other techniques; for blocking fibrinolytic activity in vitro.

Aprotinin is a multiple inhibitor of many serine proteases, including Trypsin, Plasmin, Urokinase, Kallikrein, Factor XIIa. It contains 58 amino acids. It forms stable complexes with these enzymes.

5D-62106-RUO- Aprotinin 250 mg Lyophilized Powder [Click here to view more information](#)

This product is for Research Use Only and not for use in Diagnostic procedures.

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.

Prionex® is a registered trademark of DSM IP Assets B.V."



Contact Information

Maria Waldron

Aniara Diagnostica

<http://www.aniara.com>

9054840686

Online Web 2.0 Version

You can read the online version of this press release [here](#).