

## **Now Available from Aniara Diagnostica 5D Protein Stabilizer**

*Aniara Diagnostica is pleased to announce the availability of 5D Protein Stabilizer, 10% Hydrolyzed Porcine Collagen Peptides. For Research Use Only.*

WEST CHESTER, Ohio ([PRWEB](#)) February 18, 2020 -- Do you use Prionex®? If so, you may want to try 5D's Protein Stabilizer!

Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

For Laboratory Use Only. Not for use in human or direct animal applications.

5D's Protein Stabilizer consists of hydrolyzed pig collagen peptide fraction, from dermal origin, which can be used as a protein stabilizer for lyophilized, liquid formulations or as a blocking agent for coated surfaces. 5D's Protein Stabilizer is free from cartilage, bone and plasma components. The product is sterilized for a long-term storage.

A5D-82411S (50 mL) 5D Protein Stabilizer [click here for more information](#)

A5D-82411A (100 mL) 5D Protein Stabilizer [click here for more information](#)

A5D-82411B (1000mL) 5D Protein Stabilizer [click here for more information](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.

"Prionex® is a registered trademark of DSM IP Assets B.V."



**Contact Information**

**Maria Waldron**

Aniara Diagnostica

<http://www.aniara.com>

513-342-3018

**Online Web 2.0 Version**

You can read the online version of this press release [here](#).