

## **Now available from Aniara Diagnostica 5-ELISA ADAMTS-13 Antigen, for measuring ADAMTS-13 antigen in plasma or in any biological fluid**

*Aniara Diagnostica is pleased to announce the availability of 5-ELISA ADAMTS-13 Antigen kit which is a two-site enzyme immunoassay for measuring ADAMTS-13 antigen in plasma or in any biological fluid where it must be tested.*

WEST CHESTER, Ohio ([PRWEB](#)) May 03, 2023 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

For Research Use Only. Not for Use in Diagnostic Procedures.

This sandwich ELISA is designed with polyclonal antibodies coated onto the plate for capturing ADAMTS-13 in the tested sample.

ADAMTS-13 is a metalloproteinase present in plasma at a concentration of about 500-800 ng/mL. It cleaves high molecular weight von Willebrand factor polymers to a lower molecular weight, following its release into blood circulation from endothelial cells (released from Weibel-Palade bodies).

5-ELISA ADAMTS-13 Antigen A5D-55501 [click here for more information](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.

**Contact Information**

**Maria Waldron**

Aniara Diagnostica

<http://www.aniara.com>

513-342-3018

**Online Web 2.0 Version**

You can read the online version of this press release [here](#).