

Now Available from Aniara Diagnostica, Bivalirudin Calibrator and Control

Aniara Diagnostica is pleased to announce the availability of Bivalirudin Calibrator and Control, lyophilized human plasma spiked with different concentrations of Bivalirudin, for the quality control of Bivalirudin assays.

WEST CHESTER, Ohio ([PRWEB](#)) June 26, 2019 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

This product is for Research Use Only and not for use in Diagnostic procedures.

Bivalirudin controls and calibrators are proposed for the quality control of anti-IIa clotting and chromogenic assays of Bivalirudin in plasma (HEMOCLLOT Thrombin Inhibitors and BIOPHEN DTI). These are lyophilized human plasmas, spiked with Bivalirudin at various concentrations, titrated and optimized for the assay of Bivalirudin.

Biophen Bivalirudin Control 225701-RUO 6 vials (2 sets of 3 vials) of 1 mL of human plasma supplemented at 2 different concentrations of Bivalirudin. [Click here for more information.](#)

Biophen Bivalirudin Calibrator 226701-RUO 10 vials (5 sets of 2 vials) of 1 mL of human plasma supplemented at 5 different concentrations of Bivalirudin. [Click here for more information.](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.



Contact Information

Maria Waldron

Aniara Diagnostica

<http://www.aniara.com>

5133423018

Online Web 2.0 Version

You can read the online version of this press release [here](#).